

## MAXIMUM RECOVERY DILUENT

A protective and isotonic diluent for maximum recovery of micro-organisms according to ISO 6887.

<b>Dehydrated media</b>	
Code number:	500 g: MRD20500, 5 kg: MRD25000
Colour:	Yellowish
Appearance:	Homogeneous hygroscopic powder
pH before autoclaving (25 °C):	6,8 – 7,2

**Direction:** Suspend **9,5 g** in one litre of distilled water and heat gently to dissolve the medium completely. Dispense into final containers and sterilise by autoclaving at 121 °C for 15 minutes.

<b>Prepared media</b>	
Bottled media:	100 ml: MRD30100, 500 ml: MRD30500
Tubed media:	150 x 15 mm: MRD40010 (10 ml)
Colour:	Yellowish
pH (25 °C):	6,9 – 7,1

**Direction:** Dispense the bottled media aseptically into sterile final containers. Media in tubes are ready to use.

### FORMULA in g/l

Peptones	1,0
Sodium chloride	8,5

**Note:** The typical formula can be adjusted to obtain optimal performance.

**Storage conditions:** Store the dehydrated media tightly closed in a dry place at room temperature. Store the bottled and tubed media protected from light at room temperature. Use before the expiry date on the label.

#### Quality control:

Test strains	Incubation temp: 37 °C	Growth	Incubation time: 24 h
<i>Escherichia coli</i> ATCC 25922		Good	

**References:** Straker and Stokes (1957) J. Appl. Bact. 26: 493.  
ISO 6887-1:2017

**In vitro diagnostic – for professional use only!**